SuperLine Product / Manual Catalog

SPMC-1805[Rev.0]

SuperLine **Product / Manual Catalog**



Dentium For Dentists By Dentists
Specifications are subject to change without any notice. Some products listed in this catalog are not available in the market due to pending approval.

HEAD OFFICE 501 Gyeonggi R&DB Center, 105 Gwanggyo-ro, Yeongtong-gu, Suwon-si, Gyeonggi-do, Korea [16229] T+82-31-888-5431 F+82-31-888-5430 E-mail biz_mail@dentium.com / www.dentium.com

SPMC-1805 [Rev.0]



Contents

SuperLine	SuperLine Characteristics	04
	SuperLine Fixture Specifications	05
	SuperLine Fixture	06
Surgical Components	Cover Screw	07
	GBR Healing Abutment	07
	Healing Abutment	08
	IOS Healing Abutment	10
Prosthetic Procedure 1	Dual / Combi Abutment - Abutment Level Impression	11
	Combi Abutment	12
	Dual Abutment	13
	Abutment Level Impression Components	15
	Restorative Kit	16
Prosthetic Procedure 2	Dual / Custom / Milling / Angled / Direct-Casting / Metal-Casting /	
	Temporary (Plastic & Ti) Abutment - Fixture Level Impression	17
	Fixture Level Impression Components	19
	Custom Abutment	20
	Milling Abutment	21
	Angled Abutment	23
	Direct-Casting Abutment	25
	Metal-Casting Abutment	25
	Ti-Temporary Abutment	26
	Plastic Temporary Abutment	26
Prosthetic Procedure 3	Screw Abutment - Abutment Level Impression	27
	Screw Abutment	28
	Angled Screw Abutment	29
	Screw Abutment Impression Components	30
Prosthetic Procedure 4	Overdenture Procedure - Positioner / Mini Ball / Magnetic Attachment	33
	Positioner	34
	Mini Ball Attachment	36
	Magnetic Attachment	37

Surgical Kit (Full)	39
Surgical Kit (Standard)	40
Surgical Kit (Short Implant)	41
Drill Stopper Kit	42
Surgical Instruments	43
DASK	51
Sinus Bur Kit	52
Sinus Kit	53
DASK / Sinus Bur Kit / TN Brush	54
Osteotome Kit	55
Trephine Kit	56
Instrument	57
Membrane Pins & Screws	58
Prosthetic Kit	59
Planning Kit	60
Prosthetic and Laboratory Instrument	61
Help Kit	62
Polymer Guide Kit	63
Sinus Kit	64
Ridge Expander Kit	65
Implant Guide Kit	66
Digital Full Kit	67
Digital Simple Kit	68
Surgical Manual	70
Prosthesis Manual	79

Instruments



S.L.A. Surface

S.L.A. (Sandblasting with Large grits and Acid etching) • Higher bone-to-implant contact • Faster bone formation on the surface

SuperLíne


reference: Kim H., et. al. "The Biocompatability of SLA-treated Titanium Implants" Biomed. Mater. 2008; 3(2):025011

In vivo test

SuperLine Characteristics

Joint stability & Improved strength for zirconia crown

· Abutment material: Grade 4 🔶 Alloy



· Long hex design: Improved recognition

Improved soft tissue management

Concave abutment design
 Non-coating

Improved wall thickness



Double thread & Tapered design




Double threaded tapered body design may provide good success rate in immediate loading cases.

reference: Kim et. al., "A Prospective, 1-year observational study of double-threaded tapered body dental implants with immediate loading" J Prosthet Dent 2015;114:46-51

Increased thread height and sharper



SuperLine Fixture Specifications



	Fixture (Mount Free)							
A	Platform	Diameter(Ø)	3.6	4.0	4.4	4.9	6.0	7.0
В	Body D	iameter(Ø)	3.6	4.0	4.5	5.0	5.0	5.8
c	Bevel	L: 7	1.5	1.5	1.5	1.5	1.5	1.5
C	Height (mm)	L: 8, 10, 12, 14	0	0.1	0.3	0.4	0.7	1.0
D	D Total Length(mm)				7, 8, 10	, 12, 14		
E	Abutment Interface(Ø)		3.33	3.33	3.33	3.33	3.33	3.33
F	F Thread Depth(mm)		0.38	0.40	0.45	0.50	0.55	0.60
	Cap Co	blor	Yellow	Green	Blue	Red	Orange	Violet
	Selection G	uideline	Anterior	Anterior	Premolar	Molar	Molar	Molar

SuperLine Fixture

• Cover screw is	not included in	n the pack	age				Unit: mm,	Scale 1.5 : 1
Platform	Body	L	Art. No.					
		7	FX 36 07 SW	07	08	10	12	14
		8	FX 36 08 SW					
Ø 3.6	Ø 3.6	10	FX 36 10 SW					
		12	FX 36 12 SW					
		14	FX 36 14 SW					1
		7	FX 40 07 SW	07	08	10	12	14
		8	FX 40 08 SW					
Ø 4.0	Ø 4.0	10	FX 40 10 SW	1	1			
		12	FX 40 12 SW					
		14	FX 40 14 SW	-				
		7	FX 45 07 SW	07	08	10	12	14
		8	FX 45 08 SW					
Ø 4.4	Ø 4.5	10	FX 45 10 SW	1	*			
		12	FX 45 12 SW			13		
		14	FX 45 14 SW	_				
		7	FX 50 07 SW	07	08 	10 ±0.4	12 ±0.4	14 ±0.4
		8	FX 50 08 SW					
Ø 4.9	Ø 5.0	10	FX 50 10 SW	11	11			
		12	FX 50 12 SW					
		14	FX 50 14 SW	-				1
		7	FX 60 07 SW	07	08	10	12	
	<i></i>	8	FX 60 08 SW					
Ø 6.0	Ø 5.0	10	FX 60 10 SW	13				
		12	FX 60 12 SW			13		
		7	FX 70 07 SW	07	08	10	12	
		8	FX 70 08 SW		5 T 11	° 🔁 Tu		
Ø 7.0	Ø 5.8	10	FX 70 10 SW		TE			
		12	FX 70 12 SW					

Cover Screw

Single use only

• Must sterilize prior to use

Unit: mm, Scale 1.5 : 1






% Hex driver: Use no more than 5N-cm of torque when screwing a cover screw to a fixture. If hex is worn, slot on the head of the product can be used to rotate it.

GBR Healing Abutment

Unit: mm, Scale 1.5 : 1






Diameter G/H

Healing Abutment



Diameter	0/11		AIL NO.
	2.0	2.0	HAB 40 20 20 E
Ø 4.0	2.0	3.5	HAB 40 20 35 E
Ø 4.0	3.0	5.0	HAB 40 30 50 E
	4.0	7.0	HAB 40 40 70 E
	2.0	2.0	HAB 45 20 20 E
Ø 4.5	2.0	3.5	HAB 45 20 35 E
	3.0	5.0	HAB 45 30 50 E
	4.0	7.0	HAB 45 40 70 E
	2.0	2.0	HAB 55 20 20 E
Ø 5.5	2.0	3.5	HAB 55 20 35 E
5.5 W	3.0	5.0	HAB 55 30 50 E
	4.0	7.0	HAB 55 40 70 E
	2.0	2.0	HAB 65 20 20 E
Ø 6.5	2.0	3.5	HAB 65 20 35 E
כ.ס ש	3.0	5.0	HAB 65 30 50 E
	4.0	7.0	HAB 65 40 70 E

Н



Unit: mm, Scale 1.5 : 1

Healing Abutment

Unit: mm, Scale 1.5 : 1


Diameter	G/H	Н	Art. No.
<i>0</i> ,7,5	2.0	3.5	HAB 75 20 35 E
Ø 7.5	3.0	5.0	HAB 75 30 50 E

Ø O F	2.0	3.5	HAB 85 20 35 E
Ø 8.5	3.0	5.0	HAB 85 30 50 E

Ø o r	2.0	3.5	HAB 95 20 35 E
Ø 9.5	3.0	5.0	HAB 95 30 50 E








Art. No.

IOS Healing Abutment


Diameter	G/H	Туре	Art. No.
	6.0	Hex	IHAB 40 06 H
	6.0	Non-hex	IHAB 40 06 N
Ø4.0	8.0	Hex	IHAB 40 08 H
Ø 4.0	8.0	Non-hex	IHAB 40 08 N
	10.0	Hex	IHAB 40 10 H
	10.0	Non-hex	IHAB 40 10 N


Diameter	G/H	Туре	Art. No.
	6.0	Hex	IHAB 50 06 H
	6.0	Non-hex	IHAB 50 06 N
<i>Ø Г 0</i>	8.0	Hex	IHAB 50 08 H
Ø 5.0	8.0	Non-hex	IHAB 50 08 N
	10.0	Hex	IHAB 50 10 H
	10.0	Non-hex	IHAB 50 10 N

Diameter	G/H	Туре	Art. No.
	6.0	Hex	IHAB 60 06 H
	6.0	Non-hex	IHAB 60 06 N
Ø 6.0	8.0	Hex	IHAB 60 08 H
Ø 6.0	8.0	Non-hex	IHAB 60 08 N
	10.0	Hex	IHAB 60 10 H
	10.0	Non-hex	IHAB 60 10 N



Unit: mm, Scale 1.5 : 1






Prosthetic Procedure 1


Impression Technique and Restoration Selection

Dual / Combi Abutment



* Hex driver: Use no more than 10N-cm of torque when screwing a healing abutment to a fixture. If hex is worn, slot on the head of the product can be used to rotate it.


Combi Abutment




ļ	Diameter	G/H	Туре	Art. No.	Туре	Art. No.
	1.0		CAB 45 10 SE		CAB 45 10 E	
	G AF	1.5	Short	CAB 45 15 SE	Long	CAB 45 15 E
		2.5		CAB 45 25 SE		CAB 45 25 E
Ø 4.5	3.5	SHOL	CAB 45 35 SE	Long	CAB 45 35 E	
		4.5		CAB 45 45 SE		CAB 45 45 E
		5.5		CAB 45 55 SE		CAB 45 55 E


	1.5		CAB 55 15 SE		CAB 55 15 E
	2.5		CAB 55 25 SE		CAB 55 25 E
Ø 5.5	3.5	Short	CAB 55 35 SE	Long	CAB 55 35 E
	4.5		CAB 55 45 SE		CAB 55 45 E
	5.5		CAB 55 55 SE		CAB 55 55 E

	1.5		CAB 65 15 SE		CAB 65 15 E
	2.5		CAB 65 25 SE		CAB 65 25 E
Ø 6.5	3.5	Short	CAB 65 35 SE	Long	CAB 65 35 E
	4.5		CAB 65 45 SE		CAB 65 45 E
	5.5		CAB 65 55 SE		CAB 65 55 E






5.5 + 5.5



Unit: mm, Scale 1 : 1

Dual Abutment_Hex



Diameter	G/H	Art. No.
	1.0	DAB 45 10 HE
	1.5	DAB 45 15 HE
Ø 4.5	2.5	DAB 45 25 HE
Ø 4.5	3.5	DAB 45 35 HE
	4.5	DAB 45 45 HE
	5.5	DAB 45 55 HE
		1
	1.5	DAB 55 15 HE
	2.5	DAB 55 25 HE
Ø 5.5	3.5	DAB 55 35 HE
	4.5	DAB 55 45 HE
	5.5	DAB 55 55 HE
	1.5	DAB 65 15 HE
	2.5	DAB 65 25 HE
Ø 6.5	3.5	DAB 65 35 HE
	4.5	DAB 65 45 HE
	5.5	DAB 65 55 HE

 $\begin{array}{c} 5.5 \\ 5.5 \\ 1.0 \end{array} \begin{array}{c} 5.5 \\ 1.5 \end{array} \begin{array}{c} 5.5 \\ 1.5 \end{array} \begin{array}{c} 5.5 \\ 2.5 \end{array} \begin{array}{c} 5.5 \\ 3.5 \end{array} \begin{array}{c} 5.5 \\ 4.5 \end{array} \begin{array}{c} 5.5 \\ 5.5 \end{array} \begin{array}{c} 5.5 \\ 5.5 \end{array} \end{array}$ $\begin{array}{c} 5.5 \\ 1.5 \\ 1.5 \end{array} \begin{array}{c} 5.5 \\ 2.5 \end{array} \begin{array}{c} 5.5 \\ 3.5 \end{array} \begin{array}{c} 5.5 \\ 3.5 \end{array} \begin{array}{c} 5.5 \\ 4.5 \end{array} \begin{array}{c} 5.5 \\ 5.5 \end{array} \begin{array}{c} 5.5 \\ 5.5 \end{array} \begin{array}{c} 5.5 \\ 5.5 \end{array} \end{array}$





Unit: mm, Scale 1.5 : 1

· Abutment screw is included


Dual Abutment_Non-hex








Unit: mm, Scale 1.5 : 1



Abutment Level Impression Components

Unit: mm, Scale 1 : 1

Comfort Cap I Snap on			
Туре	Diameter	Art. No.	
	Ø 4.5	CCC 45 CS	
Short	Ø 5.5	CCC 55 CS	
	Ø 6.5	CCC 65 CS	
	Ø 4.5	CCC 45 C	
Long	Ø 5.5	CCC 55 C	
	Ø 6.5	CCC 65 C	



Impression Coping

Diameter	Art. No.
Ø 4.5	CIC 45 L
Ø 5.5	CIC 55 L
Ø 6.5	CIC 65 L




Lab Analog		
Туре	Diameter	
	Ø 4.5	
Short	0 5 5	


	Ø 4.5	CAN 45 SL
Short	Ø 5.5	CAN 55 SL
	Ø 6.5	CAN 65 SL
	Ø 4.5	CAN 45 LL
Long	Ø 5.5	CAN 55 LL
	Ø 6.5	CAN 65 LL

Art. No.

Burn-out Cylinder

Туре	Diameter	Art. No.
	Ø 4.5	CBC 45 SL
Single	Ø 5.5	CBC 55 SL
	Ø 6.5	CBC 65 SL
	Ø 4.5	CBC 45 BL
Bridge	Ø 5.5	CBC 55 BL
	Ø 6.5	CBC 65 BL






* Note: It is recommended to keep the torque level at 25~30 N·cm to tighten the dual abutment with fixture.

Restorative Kit

Prosthetic Procedure 2

Impression Technique and Restoration Selection

Dual / Custom / Milling / Angled / Direct-Casting / Metal-Casting / Ti-Temporary / Plastic Temporary Abutment





Combi & Dual Abutment

Aut No	Lab. Components				
Art. No	Comfort Cap	Impression Coping	Analog	Burn-out	Cylinder
XSDAB 45 S XSDAB 45	CCC 45 CS CCC 45 C	CIC 45 L	CAN 45 SL CAN 45 LL	CBC 45 SL	CBC 45 BL
XSDAB 55 S XSDAB 55	CCC 55 CS CCC 55 C	CIC 55 L	CAN 55 SL CAN 55 LL	CBC 55 SL	CBC 55 BL
XSDAB 65 S XSDAB 65	CCC 65 CS CCC 65 C	CIC 65 L	CAN 65 SL CAN 65 LL	CBC 65 SL	CBC 65 BL

Unit: mm, Scale 1 : 1

Fixture Level Impression Components

· Impression coping screw is included with Impression coping.

Impression	Coping	Pick-up
mpression	coping	i icit up

Diameter	Size	Туре	Art. No.
	Short	Hex	DPU 40 11 HE
640	Short	Non-hex	DPU 40 11 NE
Ø 4.0	Long	Hex	DPU 40 15 HE
	Long	Non-hex	DPU 40 15 NE

	Short	Hex	DPU 45 11 HE
Ø 4.5	Short	Non-hex	DPU 45 11 NE
Ø 4.5	Long	Hex	DPU 45 15 HE
	Long	Non-hex	DPU 45 15 NE


	Short	Hex	DPU 55 11 HE
Ø F F	Short	Non-hex	DPU 55 11 NE
Ø 5.5	Long	Hex	DPU 55 15 HE
	Long	Non-hex	DPU 55 15 NE

	Short	Hex	DPU 65 11 HE
0 4 F	Short	Non-hex	DPU 65 11 NE
Ø 6.5	Long	Hex	DPU 65 15 HE
	Long	Non-hex	DPU 65 15 NE

Impression Coping Transfer Screw

Size	L	Art. No.
Short	23.0	DPS 11 E
Long	27.0	DPS 15 E

Analog		
Diameter	L	Art. No.
Ø 4.5	12.0	DANSE




Unit: mm, Scale 1 : 1

Non-hex Non-hex Hex



Hex




Fixture Level Impression Components


· Impression coping screw is included with Impression coping.

Impression Coping Transfer

Diameter	Size	Туре	Art. No.
	Short	Hex	DTF 40 11 HE
<i><i>α</i>₁₀</i>	Short	Non-hex	DTF 40 11 NE
Ø 4.0	Long	Hex	DTF 40 15 HE
	Long	Non-hex	DTF 40 15 NE

	Short	Hex	DTF 45 11 HE
Ø 4 F	Short	Non-hex	DTF 45 11 NE
Ø 4.5	Long	Hex	DTF 45 15 HE
	Long	Non-hex	DTF 45 15 NE





	Short	Hex	DTF 55 11 HE
Ø F F	Short	Non-hex	DTF 55 11 NE
Ø 5.5	Long	Hex	DTF 55 15 HE
	Long	Non-hex	DTF 55 15 NE


	Short	Hex	DTF 65 11 HE
<i>M</i> 6 F	Short	Non-hex	DTF 65 11 NE
Ø 6.5	Long	Hex	DTF 65 15 HE
	Long	Non-hex	DTF 65 15 NE


Impression Coping Transfer Screw

Size	L	Art. No.
Short	17.0	DTS 11 E
Long	21.0	DTS 15 E

Analog

linalog		
Diameter	L	Art. No.
Ø 4.5	12.0	DANSE






Hex

Non-hex

Hex

Non-hex




Ø 4.5 12.0 · Abutment screw is included.

Unit: mm, Scale 1.5 : 1

Custom Abutment



Diameter	G/H	Туре	Art. No.
	0.5	Hex	CDAB 45 05 HE
Ø 4.5	0.5	Non-hex	CDAB 45 05 NE
Ø 4.5	1.5	Hex	CDAB 45 15 HE
	1.5	Non-hex	CDAB 45 15 NE
	1.0	Hex	CDAB 55 10 HE
Ø 5.5	1.0	Non-hex	CDAB 55 10 NE
د.د ש	2.0	Hex	CDAB 55 20 HE
	2.0	Non-hex	CDAB 55 20 NE



Unit: mm, Scale 1.5 : 1

Milling Abutment

· Abutment screw is included.



Ø 4.0	Hex	MAB 40 105 HE
	Non-hex	MAB 40 105 NE

Diameter	G/H	Туре	Art. No.
<i>Q</i> 4 F	1.5	Hex	MAB 45 156 HE
Ø 4.5	1.5	Non-hex	MAB 45 156 NE

Diameter	G/H	Туре	Art. No.
	1.5	Hex	MAB 55 157 HE
Ø F F	1.5	Non-hex	MAB 55 157 NE
Ø 5.5	2.5	Hex	MAB 55 257 HE
	2.5	Non-hex	MAB 55 257 NE








· Abutment screw is included.

Milling Abutment

Diameter	G/H	Туре	Art. No.
	1.5	Hex	MAB 65 158 HE
	1.5	Non-hex	MAB 65 158 NE
a.c.	2.5	Hex	MAB 65 258 HE
Ø 6.5	2.5	Non-hex	MAB 65 258 NE
	3.5	Hex	MAB 65 358 HE
	3.5	Non-hex	MAB 65 358 NE



Unit: mm, Scale 1.5 : 1

Diameter	G/H	Туре	Art. No.
	2.5	Hex	MAB 75 259 HE
<i>M</i> 7 F	2.5	Non-hex	MAB 75 259 NE
Ø 7.5	3.5	Hex	MAB 75 359 HE
	3.5	Non-hex	MAB 75 359 NE




Angled Abutment_15°

· Abutment screw is included.



Diameter	G/H	Туре	Art. No.
	1.5	Hex	AAB 15 45 15 HE
	1.5	Non-hex	AAB 15 45 15 NE
Ø 4.5	2.5	Hex	AAB 15 45 25 HE
Ø 4.5	2.5	Non-hex	AAB 15 45 25 NE
	3.5	Hex	AAB 15 45 35 HE
	3.5	Non-hex	AAB 15 45 35 NE
	1		1



Diameter	G/H	Туре	Art. No.
	1.5	Hex	AAB 15 55 15 HE
	1.5	Non-hex	AAB 15 55 15 NE
Ø 5.5	2.5	Hex	AAB 15 55 25 HE
5.5 W	2.5	Non-hex	AAB 15 55 25 NE
	3.5	Hex	AAB 15 55 35 HE
	3.5	Non-hex	AAB 15 55 35 NE




Unit: mm, Scale 1.5 : 1

· Abutment screw is included.

Unit: mm, Scale 1.5 : 1

Unit: mm, Scale 1.5 : 1

Angled Abutment_25°



Diameter	G/H	Туре	Art. No.
	1.5	Hex	AAB 25 45 15 HE
	1.5	Non-hex	AAB 25 45 15 NE
Ø4.5	2.5	Hex	AAB 25 45 25 HE
Ø 4.5	2.5	Non-hex	AAB 25 45 25 NE
	3.5	Hex	AAB 25 45 35 HE
	3.5	Non-hex	AAB 25 45 35 NE

Diameter	G/H	Туре	Art. No.	
	1.5	Hex	AAB 25 55 15 HE	
	1.5	Non-hex	AAB 25 55 15 NE	
Ø F F	2.5	Hex	AAB 25 55 25 HE	
Ø 5.5	2.5	Non-hex	AAB 25 55 25 NE	(
	3.5	Hex	AAB 25 55 35 HE	
	3.5	Non-hex	AAB 25 55 35 NE	



Unit: mm, Scale 1.5 : 1



Direct-Casting Abutment





Metal-Casting Abutment

• Abutment screw is included.






* Note: It is recommended to keep the torque level at 25~30 N-cm to tighten the angled abutment with fixture.

Abutment screw is included.

Ti-Temporary Abutment



Plastic Temporary Abutment

• Abutment screw is included.

Diameter Ø 4.5

Diameter	G/H	Туре	Art. No.
<i>Ø</i> 4 F	2.0	Hex	RAB 45 20 PHL
Ø 4.5	2.0	Non-hex	RAB 45 20 PNL
Ø F F	2.0	Hex	RAB 55 20 PHL
Ø 5.5	2.0	Non-hex	RAB 55 20 PNL
0.45	2.0	Hex	RAB 65 20 PHL
Ø 6.5	2.0	Non-hex	RAB 65 20 PNL
Hex N	on-hex	2.0 Hex	12.0 Non-hex
Hex IN	on-nex	нех	Non-nex

Hex

Non-hex

Unit: mm, Scale 1.5 : 1


Unit: mm, Scale 1.5 : 1

Note: It is recommended to keep the torque level at 25~30 N-cm to tighten the temporary abutment with fixture. **Prosthetic Procedure 3**

Impression Technique and Restoration Selection


Screw Abutment





Unit: mm, Scale 1.5 : 1

Screw Abutment



Diameter	G/H	Art. No.
	1.0	SAB 45 10 E
	1.5	SAB 45 15 E
Ø 4.5	2.5	SAB 45 25 E
Ø 4.5	3.5	SAB 45 35 E
	4.5	SAB 45 45 E
	5.5	SAB 45 55 E

Diameter	G/H	Art. No.
	1.5	SAB 55 15 E
	2.5	SAB 55 25 E
Ø 5.5	3.5	SAB 55 35 E
	4.5	SAB 55 45 E
	5.5	SAB 55 55 E



G/H 3.5

SAB5535E


FX4510SW

Unit: mm, Scale 1.5 : 1




Angled Screw Abutment

· Abutment screw is included.






	Diameter	G/H	Туре	Art. No.
		1.0	Hex	ASA 45 10 HE
		1.0	Non-hex	ASA 45 10 NE
	<i>O</i> AF	2.0	Hex	ASA 45 20 HE
	Ø 4.5	2.0	Non-hex	ASA 45 20 NE
		3.0	Hex	ASA 45 30 HE
		3.0	Non-hex	ASA 45 30 NE



Diameter	G/H	Туре	Art. No.
	1.0	Hex	ASA 55 10 HE
	1.0	Non-hex	ASA 55 10 NE
Ø F F	2.0	Hex	ASA 55 20 HE
Ø 5.5	2.0	Non-hex	ASA 55 20 NE
	3.0	Hex	ASA 55 30 HE
	3.0	Non-hex	ASA 55 30 NE



Angled Screw Abutment Screw

Diameter	L	Art. No.
Ø 2.3	7.3	ASASC 20 23



Impression Coping Pick-up

Single

Bridge

Single

Bridge

Туре

Hex

Non-hex

Hex

Non-hex

Diameter

Ø 4.5

Ø 5.5

Screw Abutment Impression Components

Art. No.

SPU 45 SL

SPU 45 BL

SPU 55 SL

SPU 55 BL

Ø 4.5	Ø 4.5
Single (Hex)	Bridge (Non-hex)
Ø 5.5	Ø 5.5

Unit: mm, Scale 1.5 : 1

Single (Hex) Bridge (Non-hex)

Ø 4.5

Bridge (Non-hex)

Ø 5.5

5

Bridge (Non-hex)

Ø 4.5

Single (Hex)

Ø 5.5


Single (Hex)

Impression Coping Transfer

Diameter	Туре		Art. No.
<i>C</i> A F	Single	Hex	STF 45 SL
Ø 4.5	Bridge Non-hex		STF 45 BL
Ø 5.5	Single	Hex	STF 55 SL
c.c @	Bridge	Non-hex	STF 55 BL


Impression Coping Screw

Туре	L	Art. No.
Pick-up	15.0	SPS 09
Transfer	9.0	STS 09




Screw Abutment Impression Components

Unit: mm, Scale 1.5 : 1





Analog		
Diameter	L	Art. No.
Ø 4.5	12.0	SAN 45 L
Ø 5.5		SAN 55 L



Polishing Protector

Diameter	L	Art. No.
Ø 4.5	4.71	SPP 45 L
Ø 5.5		SPP 55 L




Ti-Retaining Screw

Diameter	L	Art. No.
Ø 2.3	5.0	SRS 18 T



Ti-Cylinder

Diameter	Туре		Art. No.
Ø 4.5	Single	Hex	STA 45 S
	Bridge	Non-hex	STA 45 B
Ø 5.5	Single	Hex	STA 55 S
	Bridge	Non-hex	STA 55 B




Screw Abutment Impression Components

Unit: mm, Scale 1.5 : 1

Ti-Temporary Cylinder

Diameter	Туре		Art. No.
Ø 4.5	Single	Hex	STC 45 SG
	Bridge	Non-hex	STC 45 BG
Ø 5.5	Single	Hex	STC 55 SG
	Bridge	Non-hex	STC 55 BG



Ø 5.5

Ø 4.5

12.0

Gold Cylinder

Diameter	Туре		Art. No.
Ø 4.5	Single	Hex	SGC 45 SL
Ø 4.5	Bridge	Non-hex	SGC 45 BL
0.5.5	Single	Hex	SGC 55 SL
Ø 5.5	Bridge	Non-hex	SGC 55 BL

Metal	Cylinder	l Co-Cr
-------	----------	---------

Diameter	Туре		Art. No.
Ø 4.5	Single	Hex	SGC 45 CSL
<u>ю</u> н.5	Bridge	Non-hex	SGC 45 CBL
() E E	Single	Hex	SGC 55 CSL
Ø 5.5	Bridge	Non-hex	SGC 55 CBL

Burn-out Cylinder

Diameter	Туре		Art. No.
Ø 4.5	Single Bridge	Hex Non-hex	SBC 45 SL SBC 45 BL
Ø 5.5	Single	Hex	SBC 55 SL
0.5 W	Bridge	Non-hex	SBC 55 BL



Ø 4.5 Ø 5.5 12.0 12.0

Prosthetic Procedure 4


Impression Technique and Restoration Type

Overdenture Procedure

Positioner / Mini Ball / Magnetic Attachment



Positioner



Positioner Abutment

Diameter	G/H	Art. No.
	1.0	PAB 35 10
	2.0	PAB 35 20
	3.0	PAB 35 30
Ø 3.5	4.0	PAB 35 40
	5.0	PAB 35 50
	6.0	PAB 35 60
	7.0	PAB 35 70



Unit: mm, Scale 1.5 : 1

Positioner Impression Coping

Diameter	L	Art. No.
Ø 4.5	4.5	PIC

Positioner Analog

Diameter	L	Art. No.
Ø 3.5	12.4	PAN





Positioner

Unit: mm, Scale 1.5 : 1

Positioner Socket Set

Туре	Art. No.
Tilting type ±10°	FSMHS
Non Tilting type $\pm 5^{\circ}$	FSMHSN




Positioner Metal Socket

Diameter	Н	Art. No.
Ø 5.2	2.3	FSMH



Positioner Abutment

Туре	Application		Art. No.
	Blue	1,000gf	MSHP
Tilting	Orange	500gf	MSMP
type ±10°	lvory	300gf	MSLP
	White	100gf	MSOP
Non	Blue	1,000gf	MSHPN
tilting	Orange	500gf	MSMPN
type	lvory	300gf	MSLPN
±5°	White	100gf	MSOPN



Positioner Block Out Spacer

Diameter	Н	Art. No.
Ø 6.5	0.5	FSMH



Positioner Core Tool

XPCT



* Note: It is recommended to keep the torque level at 25~30 N·cm to tighten the positioner abutment with fixture.

Mini Ball Attachment



Mini Ball Abutment

Diameter	G/H	Art. No.
Ø 3.3	0	BAB 35 00 18
Ø 3.5	1.0	BAB 35 10 18
Ø 3.5	2.0	BAB 35 20 18
Ø 3.5	3.0	BAB 35 30 18
Ø 3.5	4.0	BAB 35 40 18
Ø 3.5	5.0	BAB 35 50 18

Mini Ball Impression Coping

Diameter	L	Art. No.
Ø 4.0	5.0	GICA

Mini Ball Analog


Diameter	Art. No.
Ø 3.5	BANL

Socket Spacer

Diameter	Art. No.
Ø 4.05	GBIC 3 L
Ø 4.85	GBIC 2 L

Female Socket

Diameter	Retention Force	Art. No.
Ø 4.05	300~500gf	BPF 3
Ø 4.85	500~700gf	BPF 2



Unit: mm, Scale 1.5 : 1

Ø 4.0





Ø 4.85

F

(BFS 2)

3.3

0

(BNO 2)

Magnetic Attachment_Dome Type

Unit: mm, Scale 1.5 : 1


Diameter Ø 4.5 2.0 G/H 2.0 MGT4520D MKP4520D FX4510SW

Magnetic Assay

Application	Diameter	Н	Art. No.		Ø 4.5	Ø 5.5
MKP45D	Ø 4.5	2.0	MGT 45 20 D		2.5	2.5
MKP55D	Ø 5.5	2.0	MGT 55 20 D	Retention Force	350gf	700gf

Implant Keeper Diameter

Diameter	G/H	Art. No.
	1.0	MKP 45 10 D
	2.0	MKP 45 20 D
045	3.0	MKP 45 30 D
Ø 4.5	4.0	MKP 45 40 D
	5.0	MKP 45 50 D
	6.0	MKP 45 60 D
	1.0	MKP 55 10 D
	2.0	MKP 55 20 D
Ø F F	3.0	MKP 55 30 D
Ø 5.5	4.0	MKP 55 40 D
	5.0	MKP 55 50 D
	6.0	MKP 55 60 D



** Note: It is recommended to keep the torque level at 25~30 N·cm to tighten the mini ball abutment with fixture.

Ø 4.05

(BFS 3)

2.9

O

(BNO 1)

UXNF

Magnetic Attachment_Flat Type




Magnetic Assay

Application	Diameter	Н	Art. No.				
MKP 45	Ø 4.5	1.5	MGT 45 15		Ø 4.5	Ø 4.5 Ø 4.5	Ø 4.5 Ø 4.5 Ø 5.5
	Ø 4.5	2.0	MGT 45 20		1.5		11.5 11.5 11.5
MKP 55	Ø 5.5	1.5	MGT 55 15	Retention Force	Retention Force 400gf	Retention Force 400gf 450gf	Retention Force 400qf 450qf 700qf
IVINP 55	Ø 5.5	2.0	MGT 55 20		5	5 5	

Implant Keeper Diameter

Diameter	G/H	Art. No.
	1.0	MKP 45 10
	2.0	MKP 45 20
Ø 4.5	3.0	MKP 45 30
Ø 4.5	4.0	MKP 45 40
	5.0	MKP 45 50
	6.0	MKP 45 60
	1.0	MKP 55 10
	2.0	MKP 55 20
Ø 5.5	3.0	MKP 55 30
0.5	4.0	MKP 55 40
	5.0	MKP 55 50
	6.0	MKP 55 60




Unit: mm, Scale 1.5 : 1

Surgical Kit_Full




SuperLine Surgical Full Kit



* Note: It is recommended to keep the torque level at 25~30 N·cm to tighten the magnetic abutment with fixture.


Surgical Kit_Standard






Surgical Kit_Short Implant






Drill Stopper Kit

Unit: mm, Scale 1 : 1



	Ø 2	5.5 7
First Guide Drill	Ø2.2 <u>Ø1.7</u>	8 10121416
Diameter	L	Art. No.
	26	XLD 22 07
	29	XLD 22 29
Ø 2.2	31	XLD 22 31
Ø 2.2	33	XLD 22 33
	35	XLD 22 35
	41	XLD 22 41


Surgical Instruments






XDS









Final Drill

Diameter

Ø 2.9

Ø 3.35

Ø 3.85

Ø 4.4

Ø 5.4

Surgical Instruments

0.7 8 10 12	
Art. No.	Ø2.9
XFD 34 29	Ø 3.35

XFD 38 29

XFD 43 29

XFD 48 29

XFD 58 29 SW



Unit: mm, Scale 1 : 1


F	inal Drill		0.7 8 1012 14
	Diameter	L	Art. No.
	Ø 2.9		XFD 34 31
	Ø 3.35		XFD 38 31
	Ø 3.85	31	XFD 43 31
	Ø 4.4		XFD 48 31
	Ø 5.4		XFD 58 31

L

29

	Ø 2.9 Ø 2.5
-	Ø 3.35 Ø 2.85
	Ø 3.85 Ø 3.3
	Ø4.4 Ø3.8
	Ø 5.4 Ø 4.35

Final Drill	0.7	
Diameter	L	Art. No.
Ø 2.9		XFD 34 33
Ø 3.35	33	XFD 38 33
Ø 3.85	22	XFD 43 33
Ø 4.4		XFD 48 33



Surgical Instruments

Unit: mm, Scale 1 : 1

_


35

Final Drill Stopper			
Diameter	L	Art. No.	
Ø 2.9		XFD 34 35	
Ø 3.35	35	XFD 38 35	
Ø 3.85	22	XFD 43 35	
Ø 4.4		XFD 48 35	

Final Drill		0.7 8 1	0121416
	Diameter	L	Art. No.
	Ø 2.9		XFD 34 41
	Ø 3.35	41	XFD 38 41
	Ø 3.85		XFD 43 41
	Ø 4.4		XFD 48 41

ŀ	larvest Drill	Stopper	
	Diameter	L	Art. No.
	Ø 2.85		XFH 34 35
	Ø 3.3	35	XFH 38 35
	Ø 3.85		XFH 43 35
	Ø 4.4		XFH 48 35

Final Drill For	Short Drill	
Diameter	L	Art. No.
Ø 3.6		XFC 36 07
Ø 4.0		XFC 40 07
Ø 4.5	27	XFC 45 07
Ø 5.0	27	XFC 50 07
Ø 6.0		XFC 60 07
Ø 7.0		XFC 70 07



Ø 2.9 Ø 2.5

Ø 3.35 Ø 2.85

Ø 3.85 Ø 3.3

Ø 4.4 Ø 3.8

Ø 2.9 Ø 2.5





* Note: Drill speed 1,000rpm, 30~45 N·cm with irrigation

Surgical Instruments


Countersink		
Diameter	L	Art. No.
Ø 3.6	29	XCS 36 29 SW
Ø 4.0	29	XCS 40 29 SW
Ø 4.5	29	XCS 45 29 SW
Ø 5.0	29	XCS 50 29 SW
Ø 6.0	27	XCS 60 29 SW
Ø 7.0	27	XCS 70 29 SW



Unit: mm, Scale 1 : 1

Unit: mm, Scale 1 : 1

Stopper | For first guide drill, second guide drill Diameter Drilling Depth L Art. No. XLDST 14 14 4.6 12 6.6 XLDST 12 10 XLDST 10 8.6 Ø 4.45 08 10.6 XLDST 08 06 XLDST 06 12.6 04 14.6 XLDST 04 02 XLDST 02 16.6




Stopper | For final drill 3435, 3835

Diameter	Drilling Depth	L	Art. No.
	14	4.6	XFDST 14
	12	6.6	XFDST 12
	10	8.6	XFDST 10
Ø 5.14	08	10.6	XFDST 08
	06	12.6	XFDST 06
	04	14.6	XFDST 04
	02	16.6	XFDST 02

$4.6 \quad \boxed{6.6} \quad \boxed{8.6} \quad \boxed{10.6} \quad \boxed{12.6} \quad \boxed{14.6} \quad \boxed{14.6} \quad \boxed{16.6}$

Stopper | For final drill 4335, 4835

Diameter	Drilling Depth	L	Art. No.
	14	4.6	XFDST 14 L
	12	6.6	XFDST 12 L
	10	8.6	XFDST 10 L
Ø 5.14	08	10.6	XFDST 08 L
	06	12.6	XFDST 06 L
	04	14.6	XFDST 04 L
	02	16.6	XFDST 02 L



C	Condensing D	rill	
	Diameter	L	Art. No.
	Ø 3.6	33	XCD 36 33
	Ø 4.0	33	XCD 40 33
	Ø 4.5	33	XCD 45 33
	Ø 5.0	33	XCD 50 33
	Ø 6.0	31	XCD 60 31
	Ø 7.0	31	XCD 70 31



Ø 3.6

Ø 4.0

Ø 4.5

Ø 5.0

Ø 6.0

Round Bur

Diameter	L	Art. No.
Ø 2.0	33	XRB 20 33
Ø 3.0		XRB 30 33



* Note: 1. Countersink Drill & Round Bur speed 1,000rpm, 30~45 N·cm with irrigation
 2. Condensing Drill speed 20~60rpm, 30~45 N·cm with irrigation


Surgical Instruments

Unit: mm, Scale 1 : 1

Adapter | Hex 2.5mm

Diameter	L	Art. No.
	27	XFA 27 H
Hand-piece	30	XFA 30 H
	32	XFA 32 H
	24	XFA 24 W
Detabat	26	XFA 26 W
Ratchet	29	XFA 29 W
	32	XFA 32 W





Parallel Pin I For first guide drill, second guide drill


L	Art. No.
15.5	XPP 162212
20	XPP 162220 T
26	XPP 162226 T
	20

Path Pin


Diameter	L	Art. No.
Ø 1.6	18.6	XMFPA 2

Slot Driver

Туре	L	Art. No.
Databat	17	SDA 17 R
Ratchet	25	SDA 25 R




Ø 1.6 I 18.6



Surgical Instruments

Unit: mm, Scale 1 : 1

Туре	L	Art. No.
Hand-piece	25	XHD 25 H
	21	XHD 21 W
Ratchet	25	XHD 25 W
	27	XHD 27 W
Manual	26	XHD 26 T
Manual	30	XHD 30 T



Angled Hex Driver

Туре	L	Art. No.
Hand-piece	25	XAD 25 H
	21	XAD 21 W
Ratchet	25	XAD 25 W
	27	XAD 27 W
Manual	26	XAD 26 T
Manual	30	XAD 30 T


-

Drill Extension

Diameter	L	Art. No.
Ø 3.0	26	XDE

Driver

niver		
Туре	Diameter	Art. No.
Manual	Ø 12.0	XHDHT



Ø 3.0 [26


L	Ø 12.0
6	

Surgical Instruments

Unit: mm, Scale	1	:	1	
-----------------	---	---	---	--

Screw & Ball Abutment Adapter

Туре	L	Art. No.
Hand-piece	20	XMAA 1
Ratchet	21	XMA 21 W



Mini Ball Abutment Adapter

Туре	L	Art. No.
Ratchet	21	IPST 21 W







 Torque Wrench
 Scale 0.7 : 1

 XNTW



 Depth Gauge
 Scale 0.7 : 1

 XDGL
 XDGL

 ** Note: One side of Depth Gauge measures the osteotomy depth and the other side measures the gingival height from the top



Tissue Punch


of the implant.

Diameter	L	Art. No.
Ø 4.0	27	XTS 40



DASK






Sinus Elevation Instrument




Sinus Bur Kit

















DASK / Sinus Bur Kit

Unit:	mm,	Scale	1	:	1	
-------	-----	-------	---	---	---	--

DASK Drill

Туре	DASK Drill #	Art. No.
C	DASK Drill #1	XRT 33 20 35
Crestal Approach	DASK Drill #2	XRT 37 20 35
Approach	DASK Drill #3	XED 33 10 35 D
	DASK Drill #4	XRT 06 40 25
Lateral Approach	DASK Drill #5	XRT 08 40 25
Approach	DASK Drill #6	XST 08 30 25



% Note: Drill speed 800 to 1,200rpm, 30~45N·cm with irrigation

Ø 8.0 06.6 08.0 01.53

Stopper | For XRT332035, XRT372035, XED331035D


Diameter	Drilling Depth	L	Art. No.
	08	10.6	XFDST 08
Ø F 14	06	12.6	XFDST 06
Ø 5.14	04	14.6	XFDST 04
	02	16.6	XFDST 02



TN Brush

TN Brush

Diameter	Art. No.
Ø0.8	TN-0.8
Ø1.0	TN-1.0
Ø1.3	TN-1.3



Ø 1.0 TN-1.0

Ø1.3 TN-1.3

Osteotome Kit

Unit: mm, Scale 1 : 1

Osteotome

Osteotome compresses the bone laterally, providing denser bony interface rather than removing valuable bone from the surgical site.




Osteotome Kit

<u> </u>

0	XOFK	XOFBK
	Type A (Convex)	Type B (Concave)
	XOF 20 A1	XOF 20 B1
••E	XOF 34 A1	XOF 34 B1
-•E	XOF 38 A1	XOF 38 B1
	XOF 43 A1	XOF 43 B1
0	XOF 48 A1	XOF 48 B1

Osteotome Final drill type Scale 0.4 :1			
Туре	ØA	ØВ	Art. No.
	Ø 1.7	Ø 2.8	XOF 20 A1
	Ø 2.3	Ø 2.8	XOF 34 A1
	Ø 2.7	Ø 3.2	XOF 38 A1
Type A (Convex) (Ø 2.8	Ø 3.8	XOF 43 A1
	Ø 3.0	Ø 4.3	XOF 48 A1
	Ø 1.7	Ø 2.8	XOF 20 B1
the second second	Ø 2.3	Ø 2.8	XOF 34 B1
	Ø 2.7	Ø 3.2	XOF 38 B1
Type B (Concave)	Ø 2.8	Ø 3.8	XOF 43 B1
	Ø 3.0	Ø 4.3	XOF 48 B1




Unit: mm, Scale 1 : 1.5

Trephine Kit

Unit: mm, Scale 1 : 1

Trephine Bur

- Excellent fine cutting
- Strong engagement when attaching the trephine to cortical bone
- Cut-outs facilitates ease of harvest retrieval
- 5 scale marks on the Trephine drill from 2mm to 10mm
- Easy harvesting



Tree	م من ما م	. I∕:⊥
Ire	phine	

•			• •				1.1
1.00		-			10	•	
(9)		6	0			• 11	
Ø3.0 (2.3)	• •	94.0 (3.3)	Ø5.0 (4.2)		06.0 (5.2)	· -4	
						· []	
0		à.	0		A	· 41	
		0	0		0	- 4	
Ø7.0 (6.1)		98.0 (7.1)	Ø9.0 (8.0)		@10.0 (9.0)		
	V.			XI.			

	XIT		
	Kit Inclu	des	
XTP 24 03 XTP 64 07	XTP 34 04 XTP 74 08	XTP 44 05 XTP 84 09	XTP 54 06 XTP 94 10
XII 0407	XII 7400	XII 0405	

Adapter For screw abutment	
Туре	Art. No.
Hand-piece	XMAA1
Torque Wrench	XMA 21W

Instrument





Screw Drill

Diameter	Art. No.
Ø1.2	GMD 1228



Trephine Bur		
Outer Diameter	Inner Diameter	Art. No.
Ø 3.0	Ø 2.3	XTP 24 03
Ø 4.0	Ø 3.3	XTP 34 04
Ø 5.0	Ø 4.2	XTP 44 05
Ø 6.0	Ø 5.2	XTP 54 06
Ø 7.0	Ø 6.1	XTP 64 07
Ø 8.0	Ø 7.1	XTP 74 08
Ø 9.0	Ø 8.0	XTP 84 09
Ø 10.0	Ø 9.0	XTP 94 10



Membrane screw insertion tool		[Unit: mm, Scale 1 : 0.7]
Art. No.	_	
GST26H	× >	GST 26H
GST56		
GST85		
COD		GST 56



-



Membrane Pins

	Ø3.2 H 00.55
L	Art. No.
2.80	GMT3225

Unit: mm, Scale 1 : 2




Pin case | Pin 10ea included

L	Art. No.
125.0	GMNT



Prosthetic Kit

Unit: mm, Scale 1 : 1



XIP

21

Art. No.

XHD 15

XHD 30 T

XHD 25 W

Art. No.

XMA 21 W

IPST 21 W

XMAA 1

21

25


20



Screw

L	Art. No.
2.5	GSC 15 02
4.0	GSC 15 04
6.0	GSC 15 06
8.0	GSC 15 08
10.0	GSC 15 10

Unit: mm, Scale 1 : 2





Screw case | Screw 10ea included

Hex Driver

Туре

S/T

L/T

Torque Wrench

Туре

Torque Wrench

Mini Ball

Mount

Slot Driver

Adapter

L

15

30

25

L

21

21

20

SDA 25 R

Torque Werench Scale 0.7 : 1





Planning Kit



XPK


For Combi & Dual Abutment

Diameter	G/H	Art. No.
	1.5	PDAB 45 15
	2.5	PDAB 45 25
Ø 4.5	3.5	PDAB 45 35
	4.5	PDAB 45 45
	5.5	PDAB 45 55
	1.5	PDAB 55 15
	2.5	PDAB 55 25
Ø 5.5	3.5	PDAB 55 35
	4.5	PDAB 55 45
	5.5	PDAB 55 55
	1.5	PDAB 65 15
	2.5	PDAB 65 25
Ø 6.5	3.5	PDAB 65 35
	4.5	PDAB 65 45
	5.5	PDAB 65 55

For Angled Abutment

Angled	Diameter	G/H	Art. No.
	Ø4.5	2.0	PAAB 15 45 20
15°	94.5	4.0	PAAB 15 45 40
15	Ø 5.5	2.0	PAAB 15 55 20
	Ø 5.5	4.0	PAAB 15 55 40
	Ø4.5	2.0	PAAB 25 45 20
25°	Ø 4.5	4.0	PAAB 25 45 40
25	Ø 5.5	2.0	PAAB 25 55 20
	כ.כ ש	4.0	PAAB25 55 40

Unit: mm, Scale 1 : 1



Prosthetic and Laboratory Instrument

Unit: mm, Scale 1 : 1

Hex Driver		
Hex	L	Art. No.
1.28	13	XHD 13
	15	XHD 15
	21	XHD 21
	28	XHD 28





Reamer Guide

Application	Diameter	Art. No.
	Ø 4.5	CRG 45 L
Combi / Dual Abutment	Ø 5.5	CRG 55 L
Abutment	Ø 6.5	CRG 65 L

Application	Туре	Art. No.
Screw	Bridge	SRG BL
Abutment	Single	SRG SL

Ø 4.5	





Reamer Handle		Scale 0.5 : 1
	CRH	



Hand Wrench



Reamer

Application	Art. No.
Combi / Dual Abutment	CRM
Screw Abutment	SRM



Help Kit

Unit: mm, Scale 1 : 1



XIH

Screw Remover

L	Art. No.
25	XRFS 2 S
33	XRFS 2

Screw Tap Repair

Туре	Art. No.
Тар	XRSTR
11° Guide	XRSG 11
8 ° Guide	XRSG 8

Fixture Remover

Туре	Art. No.
-	XRFRT
11°	XRFRTF
8 °	XRFRTO

Abutment Hex Remover

L	Art. No.
20	XRHR 20
25	XRHR 25

Cover & Abutment Screw Remover

L	Art. No.
25	XRRHD

Wrench		Scale 1 : 0.7
	XRFRW	















25



Polymer Guide Kit


Unit: mm, Scale 1 : 1





PGSSK

PGSCK








XSG 34 **35** S

Cantilever I	Multi-Ready (5ea)
Т	Art. No.
2.5	XSG 34 45 C

2.5




45

0

* 1 fixed Sleeve + 3 additional Sleeve Hole

· Stone Drill		XGD 23 60	(1ea)
· Guide Pin		XGP 34 23 S	(5ea)
Additional Metal Sleeve (for Cantilever Multi-Ready)	101	XPGS 34 25 A	(5ea)

Sinus Kit



GSEK


Crestal Drill

Diameter	L	Art No.
Ø3.5	35	GRT 35 35 AS



Lateral Drill

Diameter	L	Art No.
Ø6.0	28	GRT 06 50 28
Ø6.0	40	GRT 06 50 40




Sinus Curette Scale 1:2/mm



Ridge Expander Kit


Unit: mm, Scale 1 : 1



GREK

Ridge Spreader

Diameter	L	Art No.
Ø1.3 / Ø2.5	32	GRS 13 25
Ø1.8 / Ø3.0	32	GRS 18 30
Ø2.3 / Ø3.5	32	GRS 23 35



Harvest Drill

Diameter	L	Art No.
Ø3.0	29	GHD 30 29

Disc Saw

Diameter	L	Art No.
Ø8.0	25	GDS 80 25



29

Bone Chisel

Art No.	GBC 18 45 13



130




Ø 1 1 Ø3.0

53

Unit: mm, Scale 1 : 1

Implant Guide Kit



GIGK

Guide Drill

Diameter	L	Art No.
Ø2.6	29 35	GGD 26 29 GGD 26 35

Spacer

Width	Art No.
5.0	GSP 05
7.0	GSP 07
9.0	GSP 09
11.0	GSP 11

Single Guide Pin


Width	Art No.
4.0 / 5.0	GGP 04 05
6.0 / 7.0	GGP 06 07
8.0 / 9.0	GGP 08 09
10.0 / 11.0	GGP 10 11





Unit: mm, Scale 1 : 1









Digital Full Kit

FTP35 FTP45


FFL35 FFL45

FFA28H FFA28W






Digital Simple Kit



XGSSK






SGD2306 SGD2308 SGD2310 SGD2312 SFD3606 SFD3608 SFD3610 SFD3612

Contents

SURGICAL MANUAL	Surgical Drill Sequence I	70
	Surgical Drill Sequence II	71
	Drilling Depth Guide	72
	Fixture Connection	76
	Installation Procedure & Warnings	77
	Surgical Kit Maintenance	78
PROSTHESIS MANUAL	Understanding the Implant and Prosthesis	80
	Types of Abutment	81
	Dual Abutment	82
	Combi Abutment	83
	Custom / Milling Abutment	84
	Angled Abutment	85
	Direct-Casting / Metal-Casting Abutment	86
	Temporary Abutment	87
	Screw Abutment	88
	Points to Consider in Abutment Selection	89
	Minimum Height Requirement for SuperLine Prosthetic Abutment	90
	Prosthetic Procedure 1	91
	Abutment Level - Dual Abutment	92
	Abutment Level - Combi Abutment	95
	Proshtetic Procedure 2	97
	Fixture Level - [Pick-up Type] _Dual Abutment	98
	Fixture Level - [Transfer Type]_Dual Abutment	101
	Fixture Level - [Transfer Type]_Milling Abutment	104
	Fixture Level - [Pick-up Type]_Angled Abutment	106
	Fixture Level - Direct-Casting Abutment	108
	Fixture Level - Metal-Casting Abutment	109
	Fixture Level - [Pick up Type]_Temporary Abutment	110
	Proshtetic Procedure 3	111
	Abutment Level - [Transfer Type]_Screw Abutment	112
	Cementation Repair Method (SCRP)	115
	Proshtetic Procedure 4	117
	Positioner	118
	Mini Ball Attachment	121
	Magnetic Attachment	122


Surgical Drill Sequence I

Drilling Sequence Guide (Final Drill)



During Fixture Insertion, 70N·cm Torque at 50rpm is Recommended

Countersink drill is used in cases with dense cortical bone.
If the bone density is D1~D2, it is recommended to countersink after final drill.
The actual diameter of the Countersink drill is 0.1mm larger than the fixture platform.





Determination of Fixture Top Level

Top level of fixture needs to be located 0.5mm below the marginal crestal bone level to minimize bone loss after implantation.




Depth Indication

Use the depth gauge after first drill / First guide drill to check depth of drilling.
Place the depth gauge against the wall of the osteotomy.

Surgical Drill Sequence II

Drilling Sequence Guide (Condensing Drill)




Condensing Drill

Condensing Drill speed 20~60rpm, 30~45 N·cm with irrigation
If the bone density is D1~D2, it is recommended to Condensing drill after final drill.
The actual diameter of the Condensing drill is 0.1mm larger than the fixture platform.





Early loading with precise tap drill


Tap drills for exquisite surgery to enable early loading.
Smooth and precise surgery with tap drill in type 1~2 bone
4mm of tapping depth for optimum torque control.




Drilling Depth Guide









Drilling Depth Guide





Platform: Ø7.0 / Body: Ø7.0 (1000rpm / 30~45N-cm)




Fixture Connection












By hand-piece


20rpm / 35 N·cm





Directions Using the Hand-piece / Ratchet Adapter






The Hand-piece Adapter/Ratchet Adapter must be connected firmly together with the internal hex inside the fixture

Installation Procedure & Warnings

Cover Screw



Healing Abutment





Healing Abutment connection



By Hex Driver

Healing Abutment (HAB402020L) connection in thin gingiva

Warnings

Dental Implant surgery and restoration involve complex dental procedures. Appropriate and adequate training in proper technique is strongly recommended prior to use.

- Improper medical examination and/or treatment plan can result in implant failure and/or loss of supportive bone.
- Improper initial stability and/or excessive occlusal forces during healing period may lead to osseointergration failure.
- Excessive insertion torque may lead to mechanical failure or implant biologic failure due to bone compression and necrosis.
- When forces or loads are greater than its design, implant or abutment fracture could happen. Therefore clinicians should make careful decisions with regards to clinical treatment planning to minimize the risk of fracture. Appropriate implant quantity, occlusal interface and a nightguard are essential. Potential excessive loading conditions may include the following:
- 01 Inadequate number of implants are placed.
- 02 Implant width and/or length are inappropriate for a treatment site.
- 03 Prosthesis which has excessive cantilever length due to inadequate biomechanical design
- 04 Continuous occlusal force are generated by incomplete connection between implant and abutment and/or abutment screw loosening.
- 05 Direct Casting Abutment angles are greater than 30° from the vertical axis of the implant. Direct Abutments are not for angulation.
- 06 Occlusal interferences causing excessive lateral forces
- 07 Patient parafunctions such as bruxism
- 08 Inadequate dental laboratory casting procedures
- 09 Improper prosthesis fit
- 10 Trauma from patient habits or accidents
- 11 Excessive marginal bone loss caused by inadequate bone width and/or advanced periimplantitis

Surgical Kit Maintenance

Manual Cleaning and Sterilization Procedure

It is important to use protective clothing and face shield while cleaning contaminated instruments. Always wear protective glasses, mask, gloves, etc. for your safety.

Cleaning

- 1 Rinse instruments immediately after use under running tap water (<40°C) for a minimum of one (1) minute to remove all debris including extraneous body fluids, bone debris and tissue.
- 2 Soak all instruments immediately after rinsing in an enzymatic cleaning solution* for 10 to 20 minutes (Do not soak overnight).
- * Follow manufacturer's instructions and observe recommended cleaning solution concentrations (enzymatic detergent with a pH level between 7-10 and temperature not to exceed 40°C). Do not use incompatible cleaning solutions to clean instruments.
- 3 For internal irrigation drills, use a 1mL syringe and a 25 gauge needle to clean the drill irrigation hole with a minimum of 0.2 mL of the prepared cleaning solution. Repeat this step two (2) more times for a total of three (3) rinses.
- 4 Scrub with a soft brush for a minimum of 1 (one) minute to remove any debris inside the drill irrigation hole.
- 5 Rinse the instruments under running tap water (<40°C) for a minimum of 1 minute. Use a 1mL syringe and a 25 gauge needle with a minimum of 0.2 mL of tap water to forcefully flush inside the drill irrigation hole. Repeat flushing of drill irrigation hole two (2) more times for a total of three (3) flushings.
- 6 Place instruments into an ultrasonic cleaner with neutral detergent**. Keep instruments inside the ultrasonic bath for 15 minutes using a frequency of 25-50 kHz. Ensure multiple instruments placed within the bath remain separated.
- ** Follow manufacturer's instructions and observe recommended neutral detergent solution concentrations (neutral detergent with a pH level between 7-10 and temperature not to exceed 40°C). Do not use incompatible neutral detergent solutions to clean instruments.
- 7 Rinse instruments thoroughly with running tap water (<40°C) for a minimum of 1 (one) minute until all traces of neutral detergent solution are removed. Rinse inside drill irrigation hole using a 1mL syringe and a 25 gauge needle with a minimum of 0.2 mL of tap water. Repeat rinsing drill irrigation hole two (2) more times for a total of three (3) rinses.
- 8 Gently wipe instruments with a soft lint-free cloth or place the instruments in a drying cabinet (60°C for less than 10 hours) until fully dry. Blow residual water from drill irrigation hole using a 1mL syringe and a 25 gauge needle. Visually inspect instruments in a well-lit area to ensure they are clean, dry and free of residue.
- 9 Clean instrument trays with a germicidal cleaner prior to returning instruments into Kit.
- 10 Always check for damage or corrosion after rinsing and drying.

Sterilization

Dentium recommends either the Pre-vacuum or Gravity autoclave methods for sterilization under the conditions described below. However, autoclave performance can affect the efficacy of this process. Healthcare facilities should validate their sterilization processes employing the actual equipment and operators that routinely sterilize instruments.


All autoclaves/sterilizers should be regularly validated, maintained and checked in accordance with EN 285/EN 13060, EN ISO 17665, ANSI AAMI ST79 to ensure compliance with these and related standards. Make sure packaging is suitable for steam sterilization.

Recommended Sterilization Parameters

Method-Moist Heat Sterilization	Pre-vacuum	Gravity
Set Point Temperature	132 °C	132 °C
Exposure time	4 minutes	30 minutes
Drying time	20 minutes	40 minutes

PROSTHESIS MANUAL

Understanding the Implant and Prosthesis



Biological Connection

• The tapered conical hex connection between implant and abutment interface provides hermetic sealing.


• The biological connection distributes the load to the fixture evenly. Therefore it may minimize bone loss.

• All implant diameters share the same internal connection. One abutment screw fits all abutments and fixtures.

Types of Abutment (Abutments are available in various diameters & gingival heights)

· Dual Abutment · Combi Abutment	 Abutment level
 Dual Abutment Dual Milling Abutment Angled Abutment (15°/25°) Direct-Casting Abutment Metal-Casting Abutment Temporary Abutment (Plastic & Titanium) 	Fixture level
· Screw Abutment · Angled Screw Abutment (15°/ 30°)	 Screw retained (Abutment level)
 Positioner Attachment Ball Attachment Magnetic Attachment 	 For denture use

Types of Abutment




· Straight abutments are Dual and Combi Abutment.

· Depending on the insertion angle and position of the fixture, the Angled or Direct / Metal - Casting Abutment may be used.

·The Screw Abutment can be used when prosthesis retrieval is anticipated.

Dual Abutment



· It is possible to take an impression at both fixture level and abutment level.

(A Dual Abutment may be interchanged with a Combi Abutment)

- · For abutment level impressions, the same prosthetic procedures apply to both Dual and Combi Abutments.
- · For fixture level impressions, the abutment selection takes place on the master model.
- · For fixture level impressions, a precise positioning jig for abutment may be required.
- · Either hex or non-hex abutments may be used, according to operator's preference.

* If a cement retained restoration requires retrieval, cutting a hole in the occlusal surface would allow access to the screw to permit removal.


Hex / Non-hex

	Hex	Non-hex
Positioning Jig	Unnecessary	Required
Radiograph	Unnecessary	Unnecessary

Dual Abutment Line Up (Hex / Non-hex)

Diameter	G/H	Vertical Angle
Ø4.5	1.0mm, 1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm	5°
Ø5.5	1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm	6°
Ø6.5	1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm	7°

Combi Abutment



·The Combi Abutment is used when the implant position is optimal.

· If the abutment selection is made in the mouth, gauge the thickness of mucosa with the depth gauge to measure the gingival height thus allowing the appropriate abutment height.

· The Impression is taken with the snap cap.

• When using the Combi Abutment, it remains in the mouth after the impression is taken. (Do not remove or change its position)

· Tighten abutment screw to 25 - 35 N·cm. (retighten again before seating final prosthesis).

* If the Combi Abutment is too long it can be adjusted 1.5mm to the bottom of the laser mark on the vertical stack of the abutment. The Combi Abutment has a short analog for the 1.5mm adjustment.

* A resin jig can be made to record the reduction if reduced more the 1.5mm.

Combi Abutment Line Up

Diameter	G/H	Vertical Angle
Ø4.5	1.0mm, 1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm	5°
Ø5.5	1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm	6°
Ø6.5	1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm	7°

Custom / Milling Abutment





Custom Abutment

Milling Abutment

Custom Abutment

· Impression is taken at fixture level.

 \cdot When using a non-hex abutment a precise seating jig should be used.

Diameter	G/H	Туре	EIR EIR
Ø4.5	0.5mm, 1.5mm	Hey (Nep hey	
Ø5.5	1.0mm, 2.0mm	Hex / Non-hex	9

Milling Abutment

· Impression is taken at fixture level.


· When using a non-hex abutment a precise seating jig should be used.

· Either hex or non-hex abutments may be used, according to operators preference.

* If a cement retained restoration requires retrieval, cutting a hole in the occlusal surface would allow access to the screw for removal.

Diameter	G/H	Туре	C III	
Ø4.0	1.0mm			
Ø4.5	1.5mm			
Ø5.5	1.5mm, 2.5mm	Hex / Non-hex		- (]
Ø6.5	1.5mm, 2.5mm, 3.5mm			
Ø7.5	2.5mm, 3.5mm			

Angled Abutment





• The Angled Abutment is recommended when the restoration path of insertion is unfavorable in either anterior or posterior sites.

 \cdot Retention force can be increased through milling process.

Angled Abutment Line Up

Diameter	G/H	Angle
Ø4.5	1.5mm 2.5mm 3.5mm	150/250
Ø5.5	1.5mm 2.5mm 3.5mm	15° / 25°



Direct-Casting / Metal-Casting Abutment





Direct-Casting Abutment

Metal-Casting Abutment

Temporary Abutment





Ti-Temporary Abutment

Plastic Temporary Abutment

Direct-Casting Abutment

· Excellent for either single or bridgework

 \cdot Used as an esthetic custom made abutment.

· Used when angulation is not ideal and a standard abutment cannot be used.

· Used when there is inadequate inter-arch distance and a standard abutment cannot be used.

• A fixture level impression is taken, and the soft tissue contours can be supported.

Diameter	G/H	Туре
Ø4.5	1.0mm	Hex / Non-hex

Metal-Casting Abutment

· Equivalent results for a fraction of the price

· Our highly affordable metal alloy replaces expensive gold to alleviate financial burden to all.

Diameter	G/H	Туре
Ø4.5	1.0mm	Hex / Non-hex



Temporary Abutment

· Temporary Abutments are available with titanium or plastic.

• The titanium abutment comes in both hex and non-hex with a gingival height of 1.0mm. • The plastic abutment comes in diameters (Ø4.5, 5.5, 6.5) with a gingival height of 2.0mm.

Abutment	Diameter	G/H	Туре
Ti-Temporary	Ø4.5	1.0mm	Hex / Non-hex
	Ø4.5		
Plastic Temporary	Ø5.5	2.0mm	Hex / Non-hex
	Ø6.5		

Screw Abutment





Screw Abutment

Agled Screw Abutment

If prosthesis repair is anticipated, use of a Screw Abutment retained prosthesis enables easy retrieval.

· Useful for connecting multiple units or when there is a preference for a screw retained prosthesis.

· Useful when respective long axes of implants differ. Each side tapers by 30° and this permits up to 60° divergence between two abutments.

· Useful when the prognosis of an adjacent restoration is not ideal thus permitting easy retrieval and modification of the restoration.

Ti-Retaining Screw (1.8mm - body diameter)

Can minimize screw loosening due to increased approximal space.Can endure various kinds of masticatory force.

Ø 1.8

Ø 2.3

Screw Abutment Line Up

Diameter	G/H
Ø4.5	1.0mm, 1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm
Ø5.5	1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm

Angled Screw Abutment Line Up

Diameter	G/H	Angle
Ø4.5	1.0mm 2.0mm 3.0mm	20°
Ø5.5	1.0mm 2.0mm 3.0mm	20





Points to Consider in Abutment Selection

Considerations in Selecting an Abutment

Esthetic requirement
Implant angulation
Implant location
Fixture installation depth (Gingival height)
Interarch distance
Prosthesis type
Dentist & dental technician's preference

Impression of Implant

According to the case the impression can be taken at abutment or fixture level.

Fixture Level

Dual Abutment Dual Milling Abutment Angled Abutment (15° / 25°) Direct-Casting Abutment Metal-Casting Abutment Temporary Abutment (Plastic & Titanium)


Abutment Level

- 1. Dual Abutment 2. Combi Abutment 3. Screw Abutment
- 4. Angled Screw Abutment (20°)

Abutment Impression Recommendation

Abutment	Туре	Imprssion
Dual Abutment	Cementation type, screw-cementation type	Fixture level impression or abutment level impression
Combi Abutment	Cementation type	Abutment level impression
Angled Abutment	Cementation type, screw-cementation type	Fixture level impression
Screw Abutment	Screw retained type	Abutment level impression
Direct-Casting Abutment	Cementation type, screw-cementation type	Fixture level impression
Metal-Casting Abutment	Cementation type, screw-cementation type	Fixture level impression
Dual Milling Abutment	Cementation type, screw-cementation type	Fixture level impression

Minimum Height Requirement for SuperLine Prosthetic Abutment



* Diagram above indicates the minimum height required for SuperLine prosthetic abutment.

Maximum Amount of Reduction Allotted for SuperLine

Combi Abutment

Eliminate 3.0mm from the top level Combi Abutment (laser marking:1.5mm)
 Caution _ Damage may be caused to the screw if the abutment is reduced to less than 2.5mm above the gingival height.

Dual Abutment

· Preparation of the abutment top is possible as follows.

Gingival Height	Preparable Amount
1.5mm	2.0
2.5mm	3.0
3.5mm	4.0
4.5mm	5.0
5.5mm	6.0

Angled Abutment & Milling Abutment

· Preparation of the abutment top is possible as follows.

Direct-Casting Abutment & Metal-Casting Abutment

· Required minimum abutment height: at least 5.5mm above the Fixture top.


Direct-Casting Abutment & Metal-Casting Abutment

• The Screw Abutment cannot be modified, however the Casting Abutment can be modified for interarch distance, taking reduction into consideration of the height of the retaining screw.

Prosthetic Procedure 1

Impression Technique and Restoration Selection

Dual / Combi Abutment



Abutment Level_Dual Abutment



Chairside







Dual Abutment (Hex / Non-hex)






Retighten after 15 minutes Tighten it to 25~30N⋅cm.

Impression taking



Seat the plastic cap over the abutment.





Fabrication of provisional restoration or insertion of comfort cap

Abutment Level_Dual Abutment

Clinical Procedure Gold Crown Porcelain Crown T. **Final Restoration** Lab Analog Connection Cylinder Crown Wax-up Cementation Type

LabSide




Insertion of Combi Abutment analog

into impression




Make sure analog seats securely into the

impression cap. (line up the flat side of analog to the flat side of the cap)



Soft tissue model





Fabrication of master cast

Seat burn-out cylinder securely into analog.








Fabrication of metal framework

[Multiple Units]







Abutment Level_Dual Abutment





Trimming of the extended margin by using the rubber wheel

Metal framework after removal of "Lip"



Metal framework

Access hole is made when burn-out

cylinder is used to do the wax-up.

Cemented Retained Prosthesis).

SCRP: Once an access hole has been created, it can be converted to a SCRP (Screw &



[Multiple Units]

Reamer is used to eliminate "Lip" caused by 'snap-on' mechanism.



Porcelain build-up

Final prosthesis



Trim extended margin by rubber wheel




Metal framework after removal of "Lip"



Metal framework and reamer

Metal framework

Extended margin around the metal framework due to 'snap-on' mechanism



Eliminate the lip remnant caused by 'snap-on' mechanism by reamer.




Final prosthesis

Abutment Level_Combi Abutment

[Multiple Units]

Chairside







Second stage surgery (uncovering)

Following the 2nd stage surgery, soft tissue is healed around the Healing Abutment. Healing Abutment should be selected according to the size of abutment.











Snap-on the plastic impression coping with the same sized diameter abutment



Extended margin around the metal

framework due to 'snap-on' mechanism

Impression taking


Impression taking


Inner-surface of impression

Chairside









Soft tissue model

Fabrication of master cast

Placement of burn-out cylinder



modify burn-out cylinders to its proper height.













Abutment Level_Combi Abutment





[Multiple Units]

Connect the plastic bar in the middle of the trimmed burn-out cylinders to help support the resin pattern. Wax pattern may have shrinkage.





Trimming the extended margin Metal framework and reamer with a rubber wheel



Metal Framework after removal of "Lip"




Wax-up

Metal coping adaptation (Completed framework)



Porcelain build-up final prosthesis

Chairside



* If the combi analog is trimmed due to limited inter-occlusal space in the lab, make a reduction jig. Then a slight modify of the abutment in the oral cavity may be necessary to the height of the jig.




Insertion of final prosthesis and occlusal adjustment


Prosthetic Procedure 2

Impression Technique and Restoration Selection

Dual / Custom / Milling / Angled / Direct-Casting / Metal-Casting / Ti-Temporary / Plastic Temporary Abutment



Fixture Level [Pick-up Type]_Dual Abutment



Fixture Level [Pick-up Type]_Dual Abutment

[Multiple Units]

Clinical Procedure



LabSide







Connect lab analog with impression coping.

with

Fabrication of master cast



After surveying abutment milling is

possible if necessary.



Fabrication of positioning jig







Fabrication of the cap with pattern resin

Wax-up

Metal framework

Fixture Level [Pick-up Type]_Dual Abutment





Use positioning jig to transfer the abutment in model to oral cavity then tighten it to 25~30N·cm. Retighten after 15 minutes.

[Multiple Units]



* In the process of seating the prosthesis, the prosthesis can be rebounded by gingival tissue. In this case it is advised to apply occlusal load on the prosthesis for 10~15 minutes.

SCRP-Labside





Formation of access hole with long transfer coping screw



Metal framework



Final prosthesis



SCRP- Chairside

Use positioning jig to transfer the abutment in model to oral cavity then tighten it to 25~30N·cm. Retighten after 15 minutes.




Insertion of final prosthesis and adjustment of occlusion

* In the process of seating the prosthesis, the prosthesis can be rebounded by gingival tissue. In this case it is advised to apply occlusal load on the prosthesis for 10~15 minutes.

Fixture Level [Transfer Type]_Dual Abutment

[Multiple Units]











Second stage surgery (Uncovering)

Soft tissue formed around Healing Abutment


Transfer type impression coping














Impression taking



Inner surface of the impression

Fixture Level [Transfer Type]_Dual Abutment



LabSide




Impression coping and analog connection. And insert impression coping into the impression.




Make sure the impression coping is fully seated into the impression



Fabrication of master cast





Soft tissue condition after the of impression coping



Soft tissue model

Selection of Dual Abutment of proper diameter and gingival height



Verify by surveying the selected abutment. (Milling of the abutment is possible if necessary)



Fabrication of positioning jig

Fixture Level [Transfer Type]_Dual Abutment

[Multiple Units]







Seat the cap with pattern resin

Completion of wax-up





framework with porcelain





SCRP-Labside





Make an access hole in the resin cap by using the long open tray transfer screw.

Completed wax-up

SCRP- Chairside



Final prosthesis built up on the framework with porcelain



Use positioning jig to transfer the abutment in model to oral cavity then tighten it to 25~30N·cm. Retighten after 15 minutes.



Insertion of final prosthesis and occlusal adjustment. Place wax into opening of the abutment prior to sealing with composite.

* In the process of seating the prosthesis, the prosthesis can be rebounded by gingival tissue. In this case it is advised to apply occlusal load on the prosthesis for 10~15 minutes.












Fixture Level [Transfer Type]_Milling Abutment



Chairside





Placement of Healing Abutment

Impression taking



Injecting of impression material



[Single Unit]

Placement of impression coping with the same diameter as Healing Abutment

Impression coping formation on the inside of impression is observable. (Traces of impression coping on the inner surface of impression)

Labside





Fixture Level [Transfer Type]_Milling Abutment



[Single Unit]

Impression coping and analog connection. And insert impression coping into the impression.

Soft tissue model

Master cast






Selection of appropriate Dual Milling Abutment














Metal framework

Chairside



Final prosthesis



abutment in model to oral cavity then

tighten it to 25~30N·cm. Retighten after 15 minutes.


Insertion of final prosthesis and occlusal adjustment

* In the process of seating the prosthesis, the prosthesis can be rebounded by gingival tissue. In this case it is advised to apply occlusal load on the prosthesis for 10~15 minutes.





Fixture Level [Pick-up Type]_Angled Abutment



510

St .

OX.

Injecting impression material

Removed impression

Placement of pick-up impression

Unscrew, then remove the impression.

Chairside







Impression taking (individual tray with holes)

Laboratory Procedure




Fixture Level [Pick-up Type]_Angled Abutment

[Single Unit]

Labside








Impression coping with analog connections

of master model

Unscrew then separate impression from the model.









Master cast


Select an Angled Abutment.



Modification of Angled Abutment & fabrication of positioning jig




Fabrication of pattern resin cap





Metal or zirconia framework







Insertion of final prosthesis and occlusal adjustment











Abutment using positioning jig










LabSide



Placement of Direct-Casting Abutment



Fabrication of pattern resin cap



Wax-up



Metal framework



Final prosthesis




Insertion of custom abutment using positioning jig




Insertion of final prosthesis and occlusal adjustment


Fixture Level_Metal-Casting Abutment

[Single Unit]



LabSide








Fabrication of positioning jig

Placement of Metal-Casting Abutment








Metal framework

Chairside





Final prosthesis



Insertion of custom abutment using


positioning jig



Insertion of final prosthesis and occlusal adjustment



Fabrication of pattern resin cap








Wax-up



Fixture Level [Pick-up Type]_Temporary Abutment



<Using Ti Abutment>

<Using Plastic Abutment>



Considering the opposing teeth before seating the Temporary Abutment, trim off the abutment as needed and complete the Temporary Abutment prosthesis with direct resin.




Prosthetic Procedure 3

Impression Technique and Restoration Selection


Screw Abutment

Abutment Level Impression



Screw-Retained Restoration

Abutment Level [Transfer Type]_Screw Abutment



LabSide






Select and seat an appropriate Screw Abutment with delivery holder.



Screw Abutment transfer copings (abutment level)





Impression taking




Inner-surface of impression



[Multiple Units]

15 minutes with Screw Abutment adaptor.



Injecting impression material



Placement of comfort cap on Screw Abutment

Abutment Level [Transfer Type]_Screw Abutment

[Multiple Units]

Laboratory Procedure



LabSide







Connecting impression coping with Screw Abutment analog

Position impression coping and analog assembly in the exact location of the

impression









Connect the Screw Abutment cylinder then tighten it with ti-retaining screw.



Consider the distance with opposing teeth, then trim cylinder to its appropriate height.



Connect the plastic bar in the middle of trimmed burn-out cylinders to help support the wax pattern. Wax pattern may have shrinkage.



Wax-up

Abutment Level [Transfer Type]_Screw Abutment

[Multiple Units]







Completion of metal framework

Cementation Repair Method (SCRP)

[Screw & Cement Retained Prosthesis]

In Light of Implant Prosthesis:

· A screw type restoration helps to simplify prosthesis repair, including insertion and removal of the prosthesis if necessary. · Cement type restoration tend to have a stable occlusion and may enhance the adaptability.

However the weak point is that it cannot be removed after permanent cementation.

· A Dual Abutment can be cemented or screw retained.

· Combi Abutments are cement retained and no occlusal hole is necessary.

In Case of Screw Loosening or when Prosthesis Repair is Needed







Unscrew, then remove the prosthesis from the oral cavity.







Tighten the prosthesis with 25~30N⋅cm by a screw driver. * It is recommended that the abutment screw is retightened after 15 minutes.





Final prosthesis



Completion of porcelain



Removal of lip remnant in the interior of

metal framework by using reamer

Insertion of final prosthesis and occlusal adjustment. Tighten with ti-retaining screw (25~30N·cm.).


In case of the following: screw loosing Prosthesis repair In order to unscrew, form access hole on

the occlusal surface using bur.











Fill the access hole with resin.

Cementation Repair Method (SCRP)

[Screw & Cement Retained Prosthesis]

Prosthesis Separation from Abutment due to Cement Loss







Remove the screw completely with screw driver and remove prosthesis from the patient's mouth.







Finish the repair and seat it inside the patient's mouth.

Adding to the Interproximal Contact Surface due to Prosthesis Loosening




Prosthesis loosening due to contact

loosening



Form access hole using bur







Position the prosthesis in the mouth and tighten the screw with 25~30N·cm, then fill up the access hole.





Tighten the prosthesis with $25 \sim 30$ N·cm with a screw driver.



Unscrew, then remove the cemented prosthesis with abutment in the oral cavity.



Insert the prosthesis in the oral cavity and screw it in. Afterwards, perform light curing, then polish the contact area. * It is recommended that the abutment screw is retightened after 15 minutes.




Prosthetic Procedure 4

Impression Technique and Restoration Type

Overdenture Procedure

Positioner / Mini Ball / Magnetic Attachment



Positioner


Chairside






Connect the Positioner Abutment onto the fixture.





Affix the impression coping on the Positioner Abutment.



Take Impression for the production of

individual tray.

Take the final impression with the prepared individual tray.

Positioner

Case 1

Chairside



Secure spaces for the female sockets.





Place the "block out spacer" on the Connect the metal socket onto the Posi-Positioner Abutment in the patient's tioner Abutment.



Apply a small amount of resin into the space created for the metal socket.



mouth.

Remove the white plastic socket (100gf) using the positioner tool and assemble with the regular plastic socket giving the desired retention force (300, 500 or





Connect the metal socket onto the Positioner Abutment.



After the impression material is set,

discard the individual tray.



After connecting the Positioner

Abutment and the impression coping

Image of the set final impression (with impression coping)





Positioner Analog



space required for the metal socket.

Insert the Positioner Analog into the embedded impression coping.

Fabrication of denture with

conventional method





Secure spaces for the female sockets.



Positioner Abutment in the patient's mouth.





Positioner

Case 1

Chairside





Create holes for the placement of the metal sockets.

Place the "block out spacer" on the Positioner Abutment in the intraoral.



holes and the female sockets.

Apply additional resin around the metal socket where there is a shortage

of resin.

Examine the interference between inner surface of the



Apply the resin into the holes and wait until it is completely set.

Apply resin around the metal socket.



Connect the metal socket onto the Positioner Abutment.



Remove the white plastic socket (100gf) using the Positioner tool and assemble with the regular plastic socket giving the desired retention force (300, 500 or 1000gf).



After polishing, the overdenture is completed.



Case 1







Secure spaces for the female sockets.

Connect the female sockets to the Mini Ball Abutments in the intra-oral.





Female sockets are placed

in the denture.

Connect the female sockets to the Mini Ball Abutments in the intra-oral.

Female sockets are placed

in the denture.

8



Position the denture in the oral cavity and wait until the resin is completely set.

Case 2



of female sockets.



Apply the resin into the holes and wait until it is completely set.



After polishing, the overdenture is completed.



After polishing, the overdenture is completed.

Chairside

8



Examine the interference between inner surface of the holes and the



Apply resin around the female sockets.








Magnetic Attachment

Chairside





After Healing Abutment removal

Connect implant keeper with fixture and tighten it with 25~30N⋅cm.



Position the magnetic assay on the implant keeper.



Secure spaces for the magnetic assays.



Implant keepers conne the fixtures



Examine the interference between inner divot of the denture and the magnets.

Case 1



Apply resin on the divot of the denture's inner surface.



Apply some of resin around the magnetic assays.

Position the denture into the mouth and wait until the resin is completely set.

After the resin is completely set, remove excess. After polishing, the overdenture is completed.



Magnetic assays are placed in the denture.

Magnetic Attachment





the magnets.



After setting, remove denture from

the mouth.

Chairside




Position the denture in the mouth and apply small amount of resin into the hole.



Wait until the resin is completely set.



After polishing, the overdenture is completed.





Add the resin around the magnets.

DENTIUM LONG-TERM CLINICAL DATA



15 Years




Pre-op

OVER A **DECADE** OF COMMITMENT TO THE **BEST PRODUCTS** FOR DENTISTS AND PATIENTS

Post-op

Final Prosthesis



5 Years

1 Years 2 months

6 Years 8 months

9 Years 6 months